


LATITUDTM | HIP SYSTEM
Freedom of Choice

Acetabular Cup System & Uncemented Stems



Cardiovascular

Orthopedics

Diagnostics

Endo-Surgery

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our cardiovascular, orthopedics, diagnostics and endo-surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.



Meril's 2,700+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.

Meril Orthopedics | HIP SYSTEM

Meril Orthopedic, a venture of Meril in association with Maxx Ortho Inc (www.maxxmed.com), is at the helm of developing and marketing innovative Orthopedic implants. Our joint replacement technologies and wide range of products make us a valuable partner to healthcare institutions in more than 40+ countries. At Meril, we have a guiding principle that the Physician-Patient-Product interaction is of utmost importance.



Acetabular Cup System

Acetabular cup system is a modular acetabular replacement system consisting of a range of Titanium coated modular shells designed for use with highly cross-linked polyethylene liners articulating with a range of dedicated cobalt chromium alloy modular heads.

Enhance the Press Fit

The design geometry of the shell allows press fit to occur just below the acetabular bone margin to assist retention and stability.



Hip Compatibility Table

SHELL SIZES	LINER SIZE	HEAD				
40	35	22	28			
42	37	22	28			
44	37	22	28			
46	40		28	32		
48	40		28	32		
50	44		28	32	36	
52	44		28	32	36	
54	44		28	32	36	
56	48		28	32	36	40
58	48		28	32	36	40
60	52			32	36	40
62	52			32	36	40
64	52			32	36	40
66	52			32	36	40
68	52			32	36	40
70	52			32	36	40

Uncemented Stems


The design philosophy for Cement less femoral stem of latitud™ total hip replacement system is adopted from well proven and successful stem design. Cement less femoral stem is provided in two different angles to meet with patient's anatomy and developed into two versions: standard and increased offsets. The surface of cement less femoral stem is coated with hydroxyapatite with plasma spraying method for osteo-integration purpose and to be used without cement. The Dual Taper Trapezoidal Design with vertical and horizontal grooves to promote stability and resist rotational and axial displacement.

The low modulus of elasticity of titanium reduces the mismatch with the cortical bone and reduces thigh pain



Features

- The Uncemented femoral stem is made from Ti-6Al-4V ELI alloy as per ASTM F136 Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications
- Stem is coated with Osprovit® Hydroxyapatite by plasma spraying method in compliance with ISO-13779-2/ BS ISO 13779-4 below the resection line
- Porosity of coating is maximum 10% and surface roughness value for HA coated surface of uncemented femoral stem is above 30 μm with adhesion strength above 15 MPa
- The Uncemented femoral stem is straight and tapered with a lateral chamfer to aid insertion
- It has both vertical and horizontal grooves to resist axial and torsional loading
- Uncemented femoral Stem neck is having a 12/14 taper trunnion for connection with femoral head
- It is available in 11 different sizes with Provision of 135° standard, 135° lateral, 125° Standard (coxavera) neck angle and polished distal section
- Low profile polished neck increases the head-neck ratio, reducing impingement and increasing the range of motion for the patient
- Low profile shoulder is bone conserving in the greater trochanteric region and allows for ease of use in smaller incision surgery
- Metaphyseal flare ensures maximum fixation and load transfer into the proximal femur
- Horizontal grooves in the proximal stem, supported by dense cancellous bone, provide axial stability preventing subsidence



Sizes

- Comes in 135° neck angle, (0 size to 10 size)
- High offset 135° neck angle – Used for lateralization and extreme morphologies, (0 size to 10 size)
- Coxa vara 125° neck angle – Used for Coxa vara indications and lateralisation defects, (0 size to 10 size)

Surgical Steps



- 1 The proximal femur is opened using the box chisel which is positioned laterally and posteriorly so that entry is in line with the femoral intramedullary canal.



- 2 The smallest(4-8mm tapered) intramedullary reamer, which has sharp tip, is mounted on the T-Handle, and used to open up the femoral canal. Care should be taken with this first reamer and if the bone quality is poor then the 8mm intramedullary reamer should be used first in its place, as it has a more rounded tip.



- 3 Next the 9mm intramedullary reamer is used.



- 4 The smallest rasp(Size 9) is used first to prepare the proximal femur. The small tommy bar is used to control version.



5 If the fit of this first rasp is unstable then the next size of intramedullary reamer is used followed by the corresponding size rasp. This ream/rasp technique is continued until the fit of the rasp is stable.



6 With the correct rasp in place, remove the rasp handle and if required trim the neck using the calcar cutter fitted onto the T-handle.



7 With the rasp still in place, the fin cutter is gently tapped home into the groove in the rasp with the teeth facing the greater trochanter.

Instrument Set

1. Acetabulum Cup System
2. Uncemented Stem



Corporate Headquarters, Vapi
Meril Healthcare Pvt. Ltd.
H1-H3, Meril Park
Survey No. 135/2/B & 174/2
Muktanand Marg, Chala
Vapi 396 191. Gujarat. India.
T: +91 260 3052 100
F: +91 260 3052 125
E: orthopedics@merillife.com

International Sales & Marketing
Meril Healthcare Pvt. Ltd.
612, B-Wing, Bonanaza
Sahar Plaza, Andheri East
Mumbai 400 059
Maharashtra. India.
T: +91 22 3935 0700
F: +91 22 4047 9717

India Sales & Marketing
Meril Healthcare Pvt. Ltd.
512, Midas, Sahar Plaza, J.B. Nagar
Andheri East, Mumbai 400 059
Maharashtra. India.
T: +91 22 4047 9797

Meril South America
DOC MED LTDA
1079 – Cep: 04077-003 - Moema
Sao Paulo, Brazil
T/F: +55 11 3624 5935/6

Meril Germany
Bornheimer Strasse 135 – 137
D – 53119 Bonn, Germany
T/F: +49 228 7100 4000/1

Meril Turkey
Imalat Ve Ticaret A.S.
6, Mimar Sinan Mah.,
Cavusbasi Cad, Ozde Sok
Aydin Eksi Is Merkezi Kat:1
Cekmekoy/Istanbul, Turkey
T: +90 53 2272 5172

