

 **Myval**[™]
— TRANSCATHETER HEART VALVE —
PRECISION AT HEART

Meril



**At the heart of life.
At the heart of precision.**






MyVal-1 Study Favourable Outcomes at 1-year*

LOW


Device-related Mortality
Incidence of Stroke
New Permanent Pacemaker

Myval THV: Designed for Precision in Outcomes





Myval THV has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Myval THV: Unique Crimping Outcome



Upon Crimping:


- V-shaped hinges on hexagonal frame fold, generating the dense bands on fluoroscopy
- Vertical connectors 'I' give rise to light bands
- Alternating V-folds & vertical connectors give Myval THV a unique 'Zebra Crossing' appearance on fluoroscopy for ease of positioning



Myval THV is recommended to be crimped over Navigator THV Balloon Delivery System prior to insertion within introducer sheath

Myval THV: Precise Placement Technique

Schematic of Myval THV - Ideal Landing Zone



Elimination of THV frame parallax promptly ensures visualisation of characteristic dark-light bands




Fluoroscopic images have been dramatized & and may not be in 1:1 ratio

Illustrations are artist's representation and should not be considered as engineering drawings or photographs

Myval THV: Ground Zero Deployment

- Shallow deployment of Myval THV with least engagement within LVOT is possible
- Optimal orthotopic anchorage of Myval THV with marginal LVOT foot-print without risk of THV migration
- Minimal infra-annular depth $\leq 3.50\text{mm}$ avoids conduction system interference (thus minimizing the need of new permanent pacemaker dependency)

Largest circumscribable diameter in Open Cell
(for all Myval THV Diameters 20mm to 32mm)



Myval THV infra-annular depth $\leq 3.50\text{mm}$ based on THV \emptyset

Myval THV: Detailed Sizing Guide

3D Annular area mm ²		270	280	290	300	310	314	320	330	340	350	360		
3D area-derived diameter mm		18.5	18.9	19.2	19.5	19.9	20.0	20.2	20.5	20.8	21.1	21.4		
% Annular area over/under	20 mm	16.4%	12.2%	8.3%	4.7%	1.3%	0.1%	-1.8%	-4.8%	-7.6%	-10.2%	-12.7%		
	21.5 mm	34.5%	29.7%	25.2%	21.0%	17.1%	16%	13%	10%	7%	4%	0.8%		
	23 mm	53.9%	48.4%	43.3%	38.5%	34.0%	32.3%	29.8%	25.9%	22.2%	18.7%	15.4%		
3D Annular area mm ²		370	380	390	400	410	415	420	430	440	450	460	470	
3D area-derived diameter mm		21.7	22.0	22.3	22.6	22.8	23.0	23.1	23.4	23.7	23.9	24.2	24.5	
% Annular area over/under	23 mm	12.3%	9.3%	6.5%	3.9%	1.3%	0.1%	-1.1%	-3.4%	-5.6%	-7.7%	-9.7%	-11.6%	
	24.5 mm	27.4%	24.1%	20.9%	17.9%	15.0%	13.6%	12.2%	9.6%	7.1%	4.8%	2.5%	0.3%	
	26 mm	43.5%	39.7%	36.1%	32.7%	29.5%	27.9%	26.4%	23.5%	20.7%	18.0%	15.4%	13.0%	
3D Annular area mm ²		480	490	500	510	520	530	540	550	560	570	580	590	
3D area-derived diameter mm		24.7	25.0	25.2	25.5	25.7	26.0	26.2	26.5	26.7	26.9	27.2	27.4	
% Annular area over/under	26 mm	10.6%	8.4%	6.2%	4.1%	2.1%	0.2%	-1.7%	-3.5%	-5.2%	-6.9%	-8.5%	-10.0%	
	27.5 mm	23.7%	21.2%	18.8%	16.5%	14.2%	12.1%	10.0%	8.0%	6.1%	4.2%	2.4%	0.7%	
	29 mm	37.6%	34.8%	32.1%	29.5%	27.0%	24.6%	22.3%	20.1%	17.9%	15.9%	13.9%	12.0%	
3D Annular area mm ²		600	610	620	630	640	650	660	670	680	690	700	710	
3D area-derived diameter mm		27.6	27.9	28.1	28.3	28.5	28.8	29.0	29.2	29.4	29.6	29.9	30.1	
% Annular area over/under	29 mm	10.1%	8.3%	6.5%	4.8%	3.2%	1.6%	0.1%	-1.4%	-2.9%	-4.3%	-5.6%	-7.0%	
	30.5 mm	21.8%	19.8%	17.8%	16.0%	14.2%	12.4%	10.7%	9.0%	7.4%	5.9%	4.4%	2.9%	
	32 mm	34.0%	31.8%	29.7%	27.7%	25.7%	23.7%	21.9%	20.0%	18.3%	16.6%	14.9%	13.3%	
3D Annular area mm ²		720	730	740	750	760	770	790	800	810	820	830	840	
3D area-derived diameter mm		30.3	30.5	30.7	30.9	31.1	31.3	31.7	31.9	32.1	32.3	32.5	32.7	
% Annular area over/under	32 mm	11.7%	10.2%	8.7%	7.2%	5.8%	4.4%	3.3%	1.8%	0.5%	-0.7%	-1.9%	-3.1%	-4.3%

Myval THV: Size Matrix


Myval THV Size Matrix & Technical Specifications	Area 314 mm ² 17.35 mm  20 mm	Area 363 mm ² 18.35 mm  21.5 mm	Area 415 mm ² 17.85 mm  23 mm	Area 471 mm ² 18.75 mm  24.5 mm
Perimeter	62.83 mm	67.54 mm	72.26 mm	76.97 mm
Native annulus area	270 - 330 mm ²	314 - 380 mm ²	360 - 440 mm ²	410 - 500 mm ²
Area-derived diameter	18.5 - 20.5 mm	20 - 22 mm	21.4 - 23.7 mm	22.8 - 25.2 mm
Native annulus size by TEE	16 - 19 mm	17.5 - 20.5 mm	18 - 22 mm	19.5 - 23.5 mm

All Myval THV Diameters (20 mm to 32 mm) are Compatible with 14Fr Python - Introducer Sheath

← Myval THV XL Sizes →

	Area 531 mm ²	Area 594 mm ²	Area 661 mm ²	Area 731 mm ²	Area 804 mm ²
18.85 mm					
	26 mm	27.5 mm	29 mm	30.5 mm	32 mm
	81.68 mm	86.39 mm	91.11 mm	95.82 mm	100.53 mm
	460 - 560 mm ²	510 - 630 mm ²	570 - 700 mm ²	630 - 770 mm ²	700 - 840 mm ²
	24.2 - 26.7 mm	25.5 - 28.3 mm	26.9 - 29.9 mm	28.3 - 31.3 mm	29.9 - 32.7 mm
	21 - 25 mm	22.5 - 26.5 mm	24 - 28 mm	25.5 - 29.5 mm	27 - 31 mm

Myval THV: Post Deployment Dimensions Chart



Myval THV Diameters (Ø)	20 mm	21.5 mm	23 mm
Total frame height	17.35 mm	18.35 mm	17.85 mm
Open cell height (53%)	9.20 mm	9.73 mm	9.46 mm
Closed cell height (47%)	8.15 mm	8.62 mm	8.39 mm
Infra-annular depth*	3.05 mm	3.20 mm	2.85 mm
Supra-annular height of closed cells	5.10 mm	5.42 mm	5.54 mm
Recommendation for coronary protection	10 mm	10 mm	10 mm

- A balloon occlusion test may be considered to assess the propensity for coronary occlusion
- Balloon diameter approximated to shortest axis of CT-derived annular diameter to be considered

* Infra-annular depth may depend on final landing zone position prior to valve deployment. The values mentioned are indicative range.

← Myval THV XL Sizes →

24.5 mm	26 mm	27.5 mm	29 mm	30.5 mm	32 mm
18.75 mm	18.85 mm	19.25 mm	20.35 mm	20.90 mm	21.14 mm
9.94 mm	9.99 mm	10.20 mm	10.79 mm	11.08 mm	11.21 mm
8.81 mm	8.86 mm	9.05 mm	9.56 mm	9.82 mm	9.94 mm
2.95 mm	3.05 mm	3.15 mm	3.35 mm	3.45 mm	3.55 mm
5.86 mm	5.81 mm	5.90 mm	6.21 mm	6.37 mm	6.39 mm
10 mm					

- Consider protection of coronary arteries with a DES especially if height of coronary ostium is <10 mm from the annular plane and in conjunction with sinus of valsalva dimensions i.e. height & diameters

Navigator THV Delivery System

Delivering TAVI Made Easy

- Myval THV is recommended to be crimped over Navigator THV Delivery System prior to insertion within introducer sheath
- The crimped valve with delivery system is then loaded through 14Fr Python – Introducer Sheath



- Navigator delivery system has a set of proximal and distal stoppers which ensure that valve crimping is precise and snug
- Visual confirmation of crimped valve can be ensured before entering the sheath to avoid any crimping errors/defects
- The stoppers prevent inadvertent migration of the valve & ensure there is no risk of valve dislodgement (embolization) during entry through the sheath or while negotiating the loaded delivery system across the aorta
- Myval THV direct crimping on the balloon makes TAVI delivery simple, intuitive and eliminates unwarranted procedural steps


Navigator THV Delivery System

Characteristic Balloon Expansion

- Navigator balloon with dual expansion ports at each end ensures rapid, simultaneous, controlled expansion (dog-boning) of distal and proximal ends
- This typical dog bone pattern of inflation steadies the valve during expansion phase, ensuring its precise annular position and deployment without any risk of valve migration
- Rapid balloon inflation, using an inflation device is possible with controlled palm thrust
- Rapid balloon deflation within 3-5 sec ensures procedural safety and compliance



Navigator THV Delivery System



Navigator Balloon Ø




Balloon Length



Balloon RBP




Volume of 75:25 Saline: Contrast to achieve stated balloon diameter




Navigator – THV Delivery System has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Navigator THV Delivery System: Proximal Assembly



Hi-flexion feature ensures tracking the THV delivery system via inner aortic arch curve thereby minimizing the risk of contralateral wall scraping



Caution: Always remember to fully un-flex the Navigator system while withdrawing


14Fr Python – Introducer Sheath

Compatible with all Myval THV Diameters (20 mm to 32 mm)

Sheath expands momentarily like a python swallowing its prey
 Conveniently allows passage of crimped Myval THV System

14Fr Entry Profile, Allows Atraumatic Percutaneous Access

High convenience for full retrievability of an un-deployed Myval THV System



Two separate, calibrated loading tubes ensure temporary opening of hemostatic valves in proximal port allowing smooth passage of crimped Myval THV System



Common Femoral Artery* Ø (mm)	Myval THV Ø (mm)
≥ 5.50 mm	20 mm, 21.5 mm, 23 mm, 24.5 mm
≥ 6.00 mm	26 mm, 27.5 mm, 29 mm
≥ 6.50 mm	30.5 mm, 32 mm

*CFA Ø must be MSCT derived. Excluding circumferential Ca²⁺

Python - Introducer Sheath has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Mammoth - OTW Balloon Catheter

Pre-dilatation is entirely operator's discretion and not mandatory




Mammoth Balloon Ø	14 mm	16 mm	18 mm	20 mm	23 mm	25 mm	28 mm	30 mm
Balloon length	← 40 mm →							
Balloon RBP	8 atm	← 6 atm →						
Volume of 75:25 saline: contrast to achieve stated diameter	8 ml	10 ml	13 ml	16 ml	23 ml	25 ml	34 ml	42 ml

Mammoth – OTW Balloon Catheter has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Pre-Dilatation & Balloon Sizing Rationale

- **Balloon diameter approximated to shortest axis of annular diameter measured on MSCT to be considered in order to -**
 - Ensure controlled crossing of THV across the diseased, narrowed annulus (AVA) (*prevent abrupt jumping into LV*)
 - Minimize unwarranted risk of AR due to pre-dilatation
 - Reduce propensity for iatrogenic damage to the conduction system
 - Simulate expansion behavior of the diseased native leaflets
 - Assess potential risk of ostial jailing and estimate coronary perfusion



AoV Annulus Cross-sectional View



~4000 Patients

Myval THV: Global Clinical Program

MyVal-1¹ (n=100)

Meril Initiated
Single arm
Prospective FIH study
29 Indian sites
STS ≥ 4 - ≤ 15

100 pt. 1-year follow-up presented
by Dr. Ravinder Singh Rao
at PCR Valves e-Course 2020

Myval European Study² (n=200)

Compare-TAVI³ (n=1,062)

LANDMARK Trial⁴ (n=768)

Myval-Global⁵ (n= 2,000)

Myval-China⁶ (n=125)

²Meril Initiated

Single arm

Retrospective EU Registry

20 EU sites

Real world STS

Data collection in progress

30-day f/up est.

@ PCR-LV 2021

³Investigator Initiated

RCT – 1:1
Myval : Sapien 3

Prospective study

10 Nordic + EU sites

Real world STS

Enrollment has started

1-year f/up est.

@ EuroPCR 2023

⁴Meril Initiated

RCT – 2:1:1
Myval : Sapien : Evolut THV series

Prospective study

50+ EU/ANZ sites

Real world STS

Currently enrolling

1-year f/up est.

@ EuroPCR 2023

⁵Meril Initiated

Single arm

Prospective study

100 Global sites

Real world STS

Under planning

First patient first visit est.

@ QIV 2021

⁶Meril Initiated

Single arm

Prospective study

10 China sites

Real world STS

Pre-study activities initiated

First patient first visit

est. @ QIV 2022



MyVal-1: Study Design

A prospective, multicentre, single-arm, open-label study of Myval THV in the treatment of severe symptomatic native aortic valve stenosis.

Total number of patients: 100

Device Sizes – 20, 21.5, 23, 24.5, 26, 27.5 and 29 mm

CLINICAL FOLLOW-UP



Dr. Samin Sharma - Chairman
New York, USA



Dr. Ashok Seth - Principal Investigator
New Delhi, India



Dr. Praveen Chandra - Co-ordinating PI
New Delhi, India



Dr. Ravinder Singh Rao - Co-ordinating PI
Jaipur, India



Dr. P. K. Goel - Scientific Advisor
Lucknow, India

Study Investigators: Samin Sharma, Ravinder Singh Rao, John Jose, Praveen Chandra, Pravin K. Goel, G. Sengottuvelu, Prashant Bharadwaj, C. N. Manjunath, P. C. Rath, Rajiv Chandrasekharan Nair, Rajpal Abhaichand, Ajit Mulasari, V. K. Ajith Kumar, Ajaykumar U. Mahajan, Ganesh Kumar, Jaspal Singh Arneja, Keyur Parikh, R. K. Jain, S. M. Sharma, B. B. Chanana, Jagdish Parikh, M. S. Hiremath, Rishi Sethi, Rony Mathew Kadavil, R. R. Mantri, Sanjay Mehrotra, Tarlochan Singh Kler, T. R. Murlidharan, Vijay Trehan, Ashok Seth
MyVal-1 study, one year results presented by Dr. Ravinder Singh Rao at PCR Valves e-Course 2020
MyVal-1 Study (CTRI/2016/11/007512)


Sharma SK, et.al. First-in-human evaluation of a novel balloon-expandable transcatheter heart valve in patients with severe symptomatic native aortic stenosis: the MyVal-1 study. *EuroIntervention*. 2020 Aug 28;16(5):421-9. PMID: 31566572. doi: 10.4244/EIJ-D-19-00413
Erratum for: *EuroIntervention*. 2020 Aug 28;16(5):421-429. PMID: 32855114. doi: 10.4244/EIJ-D-19-00413C

MyVal-1: Baseline Characteristics

Patient History

Average Age (years)	73.51 ± 7.49
Mean STS	5.12 ± 1.64
History of Coronary Artery Bypass Graft surgery	17%
History of Previous PCI	13%
History of Previous Aortic Valvuloplasty	1%
History of Cerebro-Vascular Events	3%
Peripheral Vascular Disease	3%

Based on the annular size,
a bioprosthetic valve of specific size was decided for each patient



- When novel intermediate sizes of the transcatheter heart valves were introduced, recruitment of only 10% of patients was remaining (of 10 patients, 5 patients were implanted with the intermediate-size study device)
- 23 mm transcatheter heart valve was implanted in 47% of the patient population- reflective of the Indian aortic annulus anatomy

MyVal-1: Cumulative Clinical Outcomes up to 1-Year Follow-up


Clinical Events [#]	Post-procedure (n=100)	1-Month Follow-Up (n=100)	6-Month Follow-Up (n=100)	12-Month Follow-Up (n=100)
Survival	98%	97%	91%	87%
All-cause Mortality	2%	3%	9%	13%
All Stroke	1%	2%	2%	2%
Disabling	0%	1%	1%	1%
Non-disabling	1%	1%	1%	1%
Acute Renal Failure	2%	2%	2%	2%
Life-threatening or disabling bleeding	1%	1%	1%	1%
Endocarditis	0%	0%	1%	1%
Myocardial Infarction	0%	0%	0%	1%
Vascular Complications	3%	3%	3%	3%
Repeat Hospitalization	NA	8%	11%	17%
New Permanent Pacemaker	2%*	2%	2%	2%

[#]One patient had RBBB pre-procedure


*Kappetein AP, Head SJ, Généreux P, et al. Updated standardized endpoint definitions for transcatheter aortic valve implantation: the Valve Academic Research Consortium-2 consensus document. Eur Heart J. 2012;33:2403-2418.

MyVal-1: Marked Improvement in Quality of Life Parameters

Six-minute walk test



Kansas City Cardiomyopathy Questionnaire Score




There was significant improvement in six-minute walk test and Kansas City Cardiomyopathy Questionnaire score from baseline to 12 months

MyVal-1: Study Conclusion

- **In 100 intermediate to high-risk patients of MyVal-1 study, the next-generation balloon-expandable transcatheter heart valve system demonstrated favourable clinical and haemodynamic outcomes at 12-month follow-up**
 - **87% survival and low incidence of all stroke (2%)**
 - **Low rate (2%) of new permanent pacemaker implantation**
 - High procedural success that can be attributed to precise orthotopic valve positioning
- In real-world global experience of 2500+ cases; Myval THV has been consistently demonstrating high procedural success and clinical performance
 - Unique hybrid honey-comb geometry for precise positioning and orthotopic deployment
 - Preserve THV geometry & respect patient's anatomy. Large size matrix : Conventional Ø 20, 23, 26, 29 mm, Intermediate Ø 21.5, 24.5, 27.5 mm & XL Ø 30.5, 32 mm
 - Direct THV crimping on Navigator balloon makes TAVI delivery simple, intuitive and eliminates unwarranted procedural steps
 - Compatibility of novel 14Fr Python Introducer Sheath for all Myval THV Øs (20-32 mm); with high convenience of full retrievability of an un-deployed Myval THV system

LANDMARK Trial – 50+ Sites EU+ANZ



Primary Endpoint – 30 Days


All cause mortality | All stroke | Life-threatening bleeding | Vascular complications
Acute kidney injury | Paravalvular leak | New permanent pacemakers

ECG/Echo Follow-up
Video Densitometry
Clinical Follow-up

Baseline | Post Procedure | 30 D | 1 Y | 3 Y | 5 Y
Post Procedure
Upto 10 years

Currently enrolling

Myval THV System and Components




Myval - THV



Navigator - THV Delivery System



Mammoth - OTW Balloon Catheter



Python - Introducer Sheath

Myval THV System and Components - Ordering Information

Myval - THV

Diameters	20.0 mm	21.5 mm	23.0 mm	24.5 mm	26.0 mm	27.5 mm	29.0 mm	30.5 mm	32.0 mm
Product code	MVL200	MVL215	MVL230	MVL245	MVL260	MVL275	MVL290	MVL305	MVL320

Navigator - THV Delivery System

Diameters	20.0 x 30 mm	21.5 x 30 mm	23.0 x 30 mm	24.5 x 30 mm	26.0 x 30 mm	27.5 x 35 mm	29.0 x 35mm	30.5 x 35 mm	32.0 x 35 mm
Product code	NVT20030	NVT21530	NVT23030	NVT24530	NVT26030	NVT27535	NVT29035	NVT30535	NVT32035

Mammoth - OTW Balloon Catheter

Diameters	14.0 x 40 mm	16.0 x 40 mm	18.0 x 40 mm	20.0 x 40 mm	23.0 x 40 mm	25.0 x 40 mm	28.0 x 40 mm	30.0 x 40 mm
Product code	MTV1440	MTV1640	MTV1840	MTV2040	MTV2340	MTV2540	MTV2840	MTV3040

Python - Introducer Sheath

Product code	PHT14
--------------	-------

Val-de-Crimp - Heart Valve Crimping Tool

Product code	VLDC
--------------	------

Myval, Navigator, Python, Mammoth & Val-de-crimp are registered trademarks of Meril Life Sciences Pvt. Ltd.

These products are intended for use by or under the direction of a trained healthcare practitioner only

Only qualified medical experts can give you information regarding your individual treatment. Prior to use, refer the instructions for use/IFU. Data on file at Meril Life Sciences Pvt. Ltd.

Illustrations are artist's representation and should not be considered as engineering drawings or photographs. Please check the regulatory approval status of Myval THV in your country

Myval, Navigator, Python, Mammoth & Val-de-crimp are not approved and not available for sale in USA. Myval, Navigator, Python, Mammoth & Val-de-crimp are not available for sale in France. Myval THV is not available for sale in Germany

Check availability of Myval THV in your country

Meril's Global Presence

Meril UK Pvt Ltd

C/O Crowe U.K LLP,
2nd Floor, 55 Ludgate Hill,
London, EC4m 7JW
United Kingdom
T: +447854524616

Meril Medical LLC.

Nauchnyi Proezd 19,
Moscow, Russia – 117 246.
Office - +7 495 772 7643

Meril GmbH.

Bornheimer Strasse 135-137,
D-53119 Bonn.
Germany
T: +49 228 7100 4000
F: +49 228 7100 4001

Meril, Inc.

2436 Emrick Boulevard,
Bethlehem, PA - 18020
USA
T: + 610 500 2080
F: + 610 317 1672

Meril South America

Doc Med LTDA Al. dos Tupiniquins,
1079 - Cep: 04077-003 - Moema,
Sao Paulo, Brazil
T: +55 11 3624 5935
F: +55 11 3624 5936

EU Representative.

Obelis S.A.Bd, General Wahis 53,
1030, Brussels, Belgium
T: +32 2 732 5954
F: +32 2 732 6003
E mail@obelis.net

Meril SA Pty Ltd

102, 104, S101 and S102,
Boulevard West Office Park,
142 Western Service Road,
Erf 813 Woodmead Extension 17
Sandton, Johannesburg – 2191
South Africa
T: +27 11 465-2049
F: +27 86 471 7941

Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş.
İçerenköy Mah. Çetinkaya Sok. Prestij Plaza No:28
Kat:4 Ataşehir, 34752
Istanbul, Turkey
T: +90 216 641 44 24
F: +90 216 641 44 25

Manufacturer

Meril Life Sciences Pvt. Ltd.

Survey No. 135/139, Bilakhia House,
Muktanand Marg, Chala,
Vapi 396 191
Gujarat, India
T: +91 260 240 8000

Meril China Co. Ltd.

Room 2301B, 23F, Lixin Plaza,
No 90, South Hubin Road,
Xiamen, China
T: 0086-592-5368505
F: 0086-592-5368519

Meril Life Sciences Pvt. Ltd.

301, A - Wing, Business Square,
Chakala, Andheri Kurla Road,
Andheri East, Mumbai - 400 093
India
T: + 91 22 39350700
E: askinfo@merillife.com



More to Life



**At the heart of life.
At the heart of precision.**

