










A Comprehensive Surgical Imaging Experience





ACCESS TO VIEW MORE...



BPL Medical[®] Technologies

MORE DETAILS IN IMAGES...

Flat Panel Detector Technology







Direct Deposition Csl

a-Si Sensor

An amorphous silicon (a-Si) image sensor for X-ray imaging consists of a two-dimensional pixelated structure. Each pixel contains a switching thin-film transistor (TFT) and a light-sensitive photodiode. Both these elements are fabricated with a-Si through photolithography process onto a large area glass substrate.

C-Ray Pro Plus uses cesium iodide doped with thallium (CsI: TI) as a scintillator. The CsI:TI is evaporated directly onto the sensor array, and is therefore in direct contact with the light sensing pixel elements. The needle-like structure of directly deposited CsI acts as light fibers, thereby preventing lateral spread of light and improving the Modulation Transfer Function (MTF).





Distortion-free & high-resolution images can be easily achieved using latest FPD technology in C-Ray Pro Plus C-Arm, which offers great visualization of details with large field of View.

MORE CONVICTION IN USAGES...




Olntuitive Interface

- Multiple Operating Mode : Cont. fluoro, Pulse fluoro, Boost fluoro, Snapshot & RAD mode
- Realtime Image enhancement features
- Realtime tube head temperature display
- Self-diagnostics

⊙ Imaging with precision

- Patient Registration (Normal & Emergency)
- Realtime noise reduction & edge enhancement
- LIH & cine loop recording
- Annotations & Measurements





• Seamless Connectivity

- DICOM 3.0 Send
- Image Export : JPEG & BMP
- Image Transfer : USB port & CD-Drive
- PACS connectivity (Optional)

MORE SAFETY IN OR...



110kHz HF Generator

produces ripple free DC high voltage based homogenous X-ray beam, which produces high quality images & reduced skin dose to the patient and OT staff.



QR-ABS (Quick Response ABS)

Quick response ABS selects optimum factors within secs that help to reduceprolong exposures



Pulsed Fuoroscopy

Parallel Collimator

Low-rate pulse fluoroscopy to reduce overall dose to the patients

Parallel shutter to collimate the beam as per ROI, which

reduces the overall dose to the patients

Laser Aimer Radiation free C-Arm positioning to minimise retakes during procedures




Removable Grid

Grid can be removed whenever it is not needed

MORE FLEXIBILITY IN OPERATION






Axial Rotation : 190° ±10°

Clearance: 805 mm

Arc Depth: 630 mm

Lateral Height: 1050 mm







Horizontal movement : 200mm



Arc Orbital :90 + 40



MORE APPLICATIONS IN OR...



Urology

General Surgery







MAJOR TECHNICAL SPECIFICATIONS	
Generator Frequency	110kHz
Power Output	3.6KW
Operating Modes	Continuous Fluoro, Boost Fluoro, Snap Mode, RAD Mode, Pulsed Fluoro
Image Receptor	Dynamic FPD
Size	21cmX21cm
Imaging Chain	1.5k X 1.5k
Image Memory	DFCC
Collimator	Parallel Shutter (Optional)
Laser Aimer	Yes
Power Supply	230 VAC ± 10% :50Hz Frequency

*Technical specification subject to change






9AM - 8PM (Monday - Saturday)

Designed to meet IEC 60601 safety standards, AERB regulations and BIS IS 7620 Part 1



A Global Medical Technology Company



ISO 13485:2016 CERTIFIED COMPANY ISO 9001:2015 CERTIFIED COMPANY

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED

Registered office: 11th KM, Bannerghatta Road, Arakere, Bangalore - 560076, India Phone +91 80 26484388/ 2648 4350 Toll Free 1800-425-2355 Website www.bplmedicaltechnologies.com CIN:U33110KA20212PTC067282 For Enquiries: sales.medical@bpl.in







Happier Living Everyday

©2024 BPL Medical Technologies Private Limited. All rights reserved. BPL Medical Technologies Private Limited reserves the right to make changes in product features, specifications, aesthetics and/or to discontinue the same at any time without notice or obligation