Instructions for use - English

1. Device description: NexGen® cobalt chromium coronary stent system includes L605 cobalt chromium stent mounted on rapid exchange balloon catheter...

Table 1: Device Component Description. Columns: Component, Dimensions (mm), Weight (g).

Table with 3 columns: Component, Dimensions (mm), Weight (g). Rows include stent length, delivery system, balloon inflation pressure, etc.

2. Indications: NexGen® is indicated for improving coronary luminal diameter in patients with a symptomatic ischemic heart disease...

3. Contraindications: Patients who are candidates for coronary bypass surgery, patients with aortic aneurysm, patients with aortic dissection...

4. Precautions: Patients with aortic dissection, patients with aortic aneurysm, patients with aortic regurgitation...

5. Warnings: Patients with aortic dissection, patients with aortic aneurysm, patients with aortic regurgitation...

6. Clinical Procedures: The following instructions should be followed to ensure proper use of the NexGen® coronary stent system...

7. Adverse Effects: Potential adverse events, which may be associated with the use of a coronary stent, include but are not limited to: abrupt closure, acute myocardial infarction...

8. Disposal: After usage, dispose the system as per accepted regulations for medical waste disposal...

9. In Vitro Information: Table - 9.1 NexGen® Stent Compliance Chart. Columns: Inflation Pressure, Stent Diameter, Stent Length.

10. Disclaimer of Warranty and Limitation of Remedy: There is no express or implied warranty, including without limitation any implied warranty of merchantability...

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Table 1: Device Component Description (Arabic). Columns: Component, Dimensions (mm), Weight (g).

Table with 3 columns: Component, Dimensions (mm), Weight (g). Rows include stent length, delivery system, balloon inflation pressure, etc.

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