
NEUROVASCULAR PRODUCTS










- High radial force and low foreshortening
- Available in 5 diameters and 3 different lengths. In addition tapered option available
- Can be deployed through a 0.014" rapid exchange as well OTW system.
- Can be used along with emboilc protection device as required.
- The device has radiopaque markers on two ends enabling proper placement of stent at desired location.



ETERNA Carotid Stent

ORDERING INFORMATION

Dia (Fr)	Length (cm)			Catheter
	40 cm	60 cm	80 cm	Length (cm)
4	√	✓	✓	
5	\checkmark	✓	✓	
6	\checkmark	✓	✓	120 cm
7	\checkmark	✓	✓	
8	\checkmark	✓	✓	

