


NOTTINGHAM ONE STEP CATHETER

Used to dilate the intramural tunnel and ureter in one step. May also be used to inject contrast medium to determine sub-mucosal tunneling or anesthetic jelly to aid in dislodging ureteral stones.

Sizes: 8/4 Fr, 10/5 Fr, 10/6 Fr, 10/7 Fr, 12/6 Fr, 12/7 Fr, 14/6 Fr, 14/7 Fr, 16/6 Fr, 16/7 Fr, 16/8 Fr, 18/6 Fr

Length: 70 cm

Model Name: BUNC



* To be used with 0.038" GW

INDOVASIVE PVT LTD

Corporate Office & Plant:

Survey No.: 48/3 & 48/7, Pashan Sus Rd, Sus Village,
Taluka Mulshi, Pune - 411021. MH, India.
Phone: 020-35000200

www.indovasive.co.in



Manufacturing Unit-01

Biorad House, #660, Hulimavu, Eshwari Industrial Estate,
Hulimavu, Bannerghatta Rd, Bangalore-560076. KA, India.
T: +91-80-2648 0683/ 6565 8007 F: +91-80-2648 0684
E: biorad@bioradmedisys.com

Manufacturing Unit-02

Plot No.18-P1, Survey No.180, Vemgal Industrial Area,
Parjenahalli Village, Vemgal Hobli, Kolar-563102 KA, India
E: biorad@bioradmedisys.com