



## NEPHROSTOMY BALLOON DILATOR

---

Used to dilate the musculofascial tract and renal capsule establishing a nephrostomy tract during a percutaneous procedure.

**Available in:** 30 Fr (ID) Amplatz Sheath

**Balloon Width:** 10 mm

**Balloon Length:** 10 cm

**Shaft Length:** 60 cm

**Shaft Size:** 7 Fr

**Model Name:**BCNB

### INDOVASIVE PVT LTD

**Corporate Office & Plant:**

Survey No.: 48/3 & 48/7, Pashan Sus Rd, Sus Village,  
Taluka Mulshi, Pune - 411021. MH, India.  
Phone: 020-35000200

[www.indovasive.co.in](http://www.indovasive.co.in)

Follow us on    

**Manufacturing Unit-01**

Biorad House, #660,Hulimavu,Eshwari Industrial Estate,  
Hulimavu,Bannerghatta Rd,Bangalore-560076. KA. India.  
T: +91-80-2648 0683/ 6565 8007 F: +91-80-2648 0684  
E: biorad@bioradmedisys.com

**Manufacturing Unit-02**

Plot No.18-P1,Survey No.180,Vemgal Industrial Area,  
Parjenahalli Village, Vemgal Hobli,Kolar-563102 KA,India  
E: biorad@bioradmedisys.com