




**Hernia Repair and Soft Tissue Augmentation with Strength,  
Softness & Sustainability**

## **FILAPROP™ MESH *Soft***

**Non Absorbable Soft Polypropylene Mesh**



### **Indications**

- ◆ Open & Laparoscopic Inguinal Hernia
- ◆ Open Ventral Hernia Repair
- ◆ Fascial Deficiencies
- ◆ Chest Wall Reconstruction

## Key Features

- 60% lighter weight than conventional polypropylene mesh.
- Affords excellent strength, durability & fast tissue integration.
- Soft & pliable mesh with unique knitting which facilitates natural process of wound healing.



## Technical Specifications

Pore Size	0.8 mm x 1.0 mm
Weight	45 gsm
Thickness	0.34 mm
Burst Strength	81.29 N/cm <sup>2</sup>
Sterilization Method	Ethylene Oxide



## Product Specifications

### Filaprop™ Mesh Soft Non-Absorbable Surgical Mesh

Sr. No.	Product Code	Brand Name	Description
1	PPMS611	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (6 CM X 11 CM)
2	PPMS715	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (7.6 CM X 15 CM)
3	PPMS1015	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (10 CM X 15 CM)
4	PPMS1215	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (12 CM X 15 CM)
5	PPMS1218	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (12 CM X 18 CM)
6	PPMS1515	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (15 CM X 15 CM)
7	PPMS1520	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (15 CM X 20 CM)
8	PPMS1530	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (15 CM X 30 CM)
9	PPMS3030	FILAPROP MESH SOFT	POLYPROPYLENE SOFT MESH (30 CM X 30 CM)