

MULTI PARAMETER MONITOR

BPL VIVID VUE 8



High Resolution Display

8.4 inch Color TFT Display with high pixel



RESP, TEMP & Thermal Recorder.(Optional)



Enhanced Connectivity Solution

Central nursing station, HL7 with HIS connectivity, USB for Screenshot & trend data transfer, DVI for external display, Nurse call function





Comprehensive Data Storage

72 Hours of Trend Data, 2000 groups of NiBP, 1000 Alarm/Arrhythmia recall - Associated waveform with all parameter data's, Full disclosure of all leads of ECG for 30min.



High performance Li-ion Battery

Low power consumption, with battery back-up time: 3.5 hrs (approx).



MULTI PARAMETER MONITOR **BPL** VIVID VUE 8



Product Specifications —

General			
Display	8.4"		
Resolution	800×600	pixels	
Display information	7 waveforms display		
Dimension	266mm×172mm×230mm		
Weight	<2.8kg		
Built-in battery		able Lithium ion battery 600mAH; ≥ 3.5 hrs of backup time	
Battery charge time	≤ 3.5 hou	rs	
ECG / Respiration Specifications			
Lead Mode	5 lead - st	andard, Optional: 3 lead	
Heart Rate Range & Accuracy	10 ~ 300 bpm Accuracy: ±1% or ≥1bpm ,whichever is greater		
Gain	2.5 mm/mV (×0.25), 5 mm/mV (×0.5), 10 mm/mV (×1), 20 mm/mV (×2)		
Sweep Speed	12.5 mm/s, 25 mm/s, 50 mm/s		
ST-segment	Detection Range -2.0mV to +2.0mV		
Respiration Method	Thoracic impedance		
RR Measurement range and Accuracy	1-150 breaths per minute Accuracy: ±2 bpm		
Delay of Apnea	10s, 15s, 20s, 25s, 30s, 35s, 40s, 45s, 50s, 55s, 60s		
NiBP Specifications			
Measurement method	Automatio	coscillometry	
Blood Pressure Range	Systolic:	ADULT 40 - 260mmHg PEDIATRIC 40 – 230mmHg NEONATE 40 – 130mmHg	
	MAP:	ADULT 26 - 220mmHg PEDIATRIC 26 - 183mmHg NEONATE 26 - 110mmHg	
	Diastolic:	ADULT 20 - 200mmHg PEDIATRIC 20 - 160mmHg NEONATE 20 - 100mmHg	
Pressure Transducer Accuracy:	± 3 mmHg between 0 mmHg - 2800 mmHg for operating conditions between 5°C - 40°C.		
Auto Measurement Interval	1min, 2min, 3min, 4min, 5min, 10min, 15min, 30min, 60min, 120min, 240min 480min, 960min.		
STAT Mode	5min		
# CE is Applicable only with Stranded Configuration.			

SPO2 Specifications		
	Masimo SpO2: Measurement range: 1 to 100%; accuracy: ±2% (adult/child, in non-motion state), ±3% (adult/child, in motion state) or ±3% (neonate, in motion or nonmotion state) within the measurement range of 70%-100%. 1%-69% is	
Measurement Range & Accuracy	not defined. PR measurement: 25-240 pulses/min.	
	Nellcor Sp02: Measurement range: 1 to 100%; accuracy: ±3% (adult/child, in non-motion state) or ±3% (neonate, in non-motion state) within the measurement range of 70%-100%. 0%-69% is not defined. PR measurement: 20-300 pulses/min	
Display resolution	1%	
Perfusion index (PI)	Masimo SpO2: 0.02%~20%; accuracy: not defined; Resolution: 0.01% (within 0.02%~9.99% range) BPL OxySat+ SpO2: 0%~20%; accuracy: not defined; Resolution: 0.01%	
Temperature Specifications		
Measurement range	10.0°C~50.0°C (50°F~122°F)	
Accuracy	± 0.1°C (± 0.2°F)	
Number of channels	2	
Optional - Recording Specifications		
Method	Thermal dot array	
Number of waveforms	3	
Recording paper width	50mm	
Paper length	15 m	
Paper speed	12.5 mm/s, 25 mm/s, 50 mm/s	
Recording way	Real-time recording, alarm recording	
Environmental Specification		
Operating temperature	+5°C to +40°C	
Operating humidity	15% to 85% (non-condensing)	
Operating atmospheric pressure	700hPa to 1060hPa	
Transportation and storage temperature	-20°C to +55°C	
Transportation and storage humidity	10% to 93% (non-condensing)	
Transportation and storage atmospheric pressure	500hPa to 1060hPa	
Power Specification		
Input voltage	100V-240V AC	
Frequency	50Hz/60Hz	
Earth leakage current	< 0.3 mA	
Input current	0.6A-1A	
Standard requirement	According to IEC 60601-1 and IEC 60601-1-2	

Technical specification subject to change

ISO 13485:2016 CERTIFIED COMPANY ISO 9001:2015 CERTIFIED COMPANY

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED

Regd. Office: 11th KM, Bannerghatta Road, Arakere, Bangalore - 560076, India.

Toll Free: 1800-4252355

Website: www.bplmedicaltechnologies.com For Enquiries: sales.medical@bpl.in CIN: U33110KA2012PTC067282





