

Disposable Bladeless Trocar

An Access System that enables hassle-free laparoscopic access and allows you to focus solely on your patient With a simple twisting motion, the bladeless tip separates the tissue, rather than cutting it.

The result: a smaller wound, lesser trauma and easier closure.

• Blunt Bladeless tip:


- Low Insertion force with simple twisting motion
- Separates rather than cuts along tissue fibers, pushing tissue and vessels away
- Contributes to faster tissue healing
- Makes the first entry safe

• Durable Primary Seal:

- Able to withstand pressure along with a long lasting performance
- Accommodates a wide range of instrument sizes

Internal Secondary Seal:

 Facilitates small specimen retrieval, Hernia mesh introduction and Laparoscopic suturing



Code	Product Description	Type Sleeve
ELR05	Bladeless Trocar, 5mm diameter, 100mm Length	Stability
ELR10	Bladeless Trocar, 10mm diameter, 100mm Length	Stability
ELR12	Bladeless Trocar, 12mm diameter, 100mm Length	Stability