

Meril

Healthcare

## MERIL TRAUMA NAILS

Crafted For Utmost Precision

With the goal of providing world class trauma implants of highest quality, Meril Trauma launches complete Range of Intermedullary Nailing System



Meril Healthcare Pvt. Ltd.  
H1-H3, Meril Park,  
Survey No. 135/2/B & 174/2  
Muktanand Marg, Chela  
Vapi 396 191 Gujarat, India  
T: +91 260 3052 100 F: +91 260 3052 125  
E: orthopedics@merilife.com

International Sales & Marketing  
Meril Healthcare Pvt. Ltd.  
612, B-Wing, Bonanaza  
Sahar Plaza, Andheri East  
Mumbai 400 059, Maharashtra, India  
T: +91 22 3935 0700 F: +91 22 4047 9717

India Sales & Marketing  
Meril Healthcare Pvt. Ltd.  
512, Midas, Sahar Plaza, J.B. Nagar  
Andheri East, Mumbai - 400 059,  
Maharashtra, India  
T: +91 22 4047 9767

Meril South America  
DOC MED LTDA  
1079 - Cep: 04077-003 - Moema  
Sao Paulo, Brazil  
T/F: +55 11 3624 5935/6

Meril Tibbi Cihazlar  
İmalat Ve Ticaret A.Ş. 6. Minhar Sınan Mah.  
Cavosbasi Cad. Cede Sok.Aydn Eksi Is Merkezi Kat:1  
Cekmekoy/Istanbul, Turkey  
T: +90 53 2272 5172

Meril GmbH  
Bornheimer Strasse 135 -137  
D - 53119 Bonn, Germany  
T/F: +49 228 7100 4000/1

[www.merilife.com](http://www.merilife.com)

## CLAVO™

### Titanium Intramedullary Nailing System

#### 1. Proximal Femoral Rotational Stability Nail - PFRN



#### 2. Cannulated Titanium Tibial Nail - CTTN



#### 3. Elastic Titanium Nail - ETN



#### 4. Shaft Femur Titanium Nail, Cannulated - SFTN



#### 5. Shaft Tibia Titanium Nail, Cannulated - STTN



### Features and Benefits

1. PFRN system comes with Anti-Rotation Screws to provide strength and angular stability
2. Design of the nails leads to high compression forces which is a desirable stimulus for fracture surface
3. M/L angle of 5 degree allow easy insertion of the nails
4. Bendable tips reduces stress on bone

## KET™

### Stainless Steel Intramedullary Nailing System

#### 1. SS Trochanteric Femoral Nail- SSTN



#### 2. SS Proximal Fixation Femoral Nail-SSFN



#### 3. SS Cannulated Femur Nail-SCFN



#### 4. SS Shaft Femur Nail, Cannulated- SSFN



#### 5. SS Shaft Tibia Nail, Cannulated-SSTN



#### 6. SS Cannulated Tibial Nail- SCFN



#### 7. SS Humeral Nail System Unreamed-SNHU



### Features and Benefits

1. Use of ASTM F138/139 Stainless Steel Grade steel leads to product of high grade with minimal amount of carbon and other impurities.
2. Anatomical design for best fit with long bones.
3. Design of locking screws provide greater biomechanical strength
4. Instrumentation sets in accordance to highest standards of quality for better precision